

IAB Onco

The IAB Onco test offers high-quality and affordable results generated using **TruSight Oncology 500 Illumina**. The data analysis and data interpretation are powered by **PierianDx** software and expert bioinformatics services.

Introduction

The IAB Onco represents a powerful assay enabling the comprehensive genomic profiling of DNA variants (523 genes) and RNA (55 genes) implicated in various tumor types. Complex assay (SNVs, InDels, fusions, splice variants, TMB, MSI) corresponding to a given tumor type can be evaluated in parallel from the same sample.

IAB Expertise

The IAB team's experiences with the TS0500 protocol are based on a number of optimisation experiments and precise laboratory work (**Fig. 1**). Simultaneous preparation of a good quality DNA and RNA highly is a comprehensive analysis and tumor characterisation. Using TS0500 ctDNA, we are able to detect tumor DNA in blood, even from very low DNA concentrations.

Material and Methods

Sample requirement

- Human DNA/RNA solid or tissue (including FFPE samples); blood
- RNA/DNA concentration: ≥ 5 ng/ μ l
- Minimum sample volume: 20 μ l

Kit for library preparation

- **TruSight Oncology 500 (Illumina)**
- Relevant DNA/RNA variants in various tumor types (**Fig. 2**)
- Targeting 523 genes (55 RNA genes)
- 1.94 Mb panel including MSI and TMB biomarkers
- Analytical specificity >99,99%



Figure 1 Workflow for comprehensive genomic profiling (CGP)




Figure 2 TS0500 library preparation and hybridization workflow

Sequencing

- NextSeq 500, High Output / NovaSeq X Plus for TS0500 ctDNA ([Illumina](#))
- Pair-end sequencing > 2x100 bp
- Internal standard – PhiX Spike-in
- BaseSpace Sequencing Hub for data storage and sharing

Data QC

- Guaranteed Phred Q30 > 80 % bases (sequencing mode 2x100bp)
- IAB recommends a minimum of 80 M/DNA and 20 M/RNA reads per sample (800 M reads for ctDNA)



Bioinformatic analysis

- Quality control and data filtering
- DNA-seq Analysis Somatic Variants
- Tumor Burden Mutation and MSI
- RNA-Seq Gene fusions
- PierianDX – Complete clinical genomic solution

Results

Sequencing results show excellent uniformity of NGS data output. Rapid data analysis, including interpretation, is performed using **PierianDx** software by expert bioinformaticians.

Conclusions

Complex NGS assay targeting 523 genes identifying all relevant DNA and RNA variants implicating in various solid tumor types (**Fig. 4**). The test also includes key current immuno-oncology biomarkers: microsatellite instability (MSI) and tumor mutation burden (TMB). Hybrid-capture chemistry combined with sophisticated bioinformatics leads to high analytical specificity and sensitivity.




Figure 3 Clinical report generated from Pierian DX software (Illumina)




Figure 4

About IAB

The experienced team follows IAB Onco workflow with strict quality procedures to perform the essential QC steps on your samples and generate reliable and repeatable results.

Certification

- ČSN EN ISO 9001:2016
- ČSN EN ISO 13485 ed. 2:2016

Illumina

Illumina is a global leader in sequencing and array technologies that are fueling groundbreaking advancements in life science research, translational and consumer genomics, and molecular diagnostics.